

Test report - PMFT 1000 - Breath resistance / Valve test



Operator	Date	Time
AnK	05.02.2021	15:55:01

Additional information

P&S Handels GmbH Auf der Gemarkung 1, 40625 Düsseldorf samir@ps-handel.com K-12806

Basic information

Operation parameter		Test conditions	
Script	Breath Resistance	Temperature	23.10 °C
Aerosol	Pressurized air	rel. humidity	2.96 %

Mask information

Mask ID	0145_20210205_155501
Mask los	Safetouch 1
Mask type	FFP2 / x / x

Result

Description	Main air	Differential pressure	Upper limit	Result
Inhalation 1	30.32 l/min	36.38 Pa	70.00 Pa	✔ pass
Inhalation 2	95.38 l/min	114.41 Pa	240.00 Pa	✔ pass
Exhalation	160.05 l/min	197.86 Pa	300.00 Pa	✔ pass



Calibration valid until 25.08.2021